



confidence
can be engineered

thermage®
CPT SYSTEM

confidence is a beautiful thing



BEFORE THERMAGE
Treatment by Ivan Rosales, MD



3 MONTHS POST TREATMENT




BEFORE THERMAGE
Treatment by Ivan Rosales, MD



3 MONTHS POST TREATMENT



BEFORE THERMAGE
Treatment by Bill Johnson, MD



6 MONTHS POST TREATMENT



BEFORE THERMAGE
Treatment by Solta Medical Aesthetic Center



4 MONTHS POST TREATMENT

practice confidence



Introducing the new Thermage CPT System.
Remarkable results, in greater comfort.

Our unique, industry leading technology has been the gold standard in skin tightening for years. That's why we're proud to introduce the new Thermage CPT™ System, a leap forward in Thermage treatment efficacy and comfort.

At the core of the system is technology that enables you to deliver more therapeutic heat, more comfortably – for better results.

The Thermage CPT System leverages our proprietary Comfort Pulse Technology™ and a new, vibrating handpiece:

- **Beautiful results** consistently, when you deliver more uniform heat to the treatment area
- **More comfort** during the procedure thanks to enhanced, proprietary Comfort Pulse Technology
- **Greater patient satisfaction** with a more comfortable treatment experience and improved results


“The Thermage CPT System delivers a significant improvement in comfort. This is very important, because when I perform a treatment, I want my patients to return for future procedures. With this new technology, 100% of the patients answered ‘yes’ for a future procedure.”

Dr. Hector Leal, Monterrey, Mexico


comfort

Significantly more comfort for patients

Our new energy delivery algorithm interweaves pulsed radiofrequency delivery with cooling bursts mimicking *Transcutaneous Electrical Nerve Stimulation* (“TENS”) therapy, to disrupt the brain’s neural response and improve comfort during treatment.



Improved Comfort



■ Current Technology
■ Thermage CPT System *Based on clinical study data

A new vibrating handpiece offers a clinically proven improvement in patient comfort.

In controlled, split-face clinical studies, comparing treatment of over 40 patients with vibration to treatment without vibration, all patients preferred being treated with vibration.


Additionally, significant improvements in average comfort scores were achieved at equivalent maximum tissue temperatures.

Our vibration technology leverages the **Gate Theory of Pain Mitigation** (Melzack and Wall, 1965). This theory states that applying “touch” to an area of discomfort can stimulate large diameter nerve fibers. These large diameter fibers then block small diameter fibers from carrying discomfort impulses to the brain.

efficacy

Our redesigned treatment tips deliver enhanced heat distribution and uniformity, which lead to improved volumetric heating. A greater volume of tissue is heated to target temperatures and the heating is more evenly distributed across the treatment area. Data from our clinical study demonstrates that the Thermage CPT System heats **four times more** tissue to therapeutic temperatures than current technologies.

Improved Heat Distribution & Uniformity



A comparison of pulses on *in-vivo* abdomen tissue demonstrates the superiority of the Thermage CPT System in delivering uniform, deep heating to the tissue

patient satisfaction

Patients will be happier with the results and the treatment experience with the Thermage CPT System.

100% of study patients preferred treatment with the Thermage CPT System vs. treatment with prior technologies.

Average patient satisfaction scores 3 months post treatment were 4.2 out of 5 (1=dissatisfied, 5=very satisfied).



thermage®

All Before and After photos are un-retouched and represent patients who have had significant results. Results vary depending on age and skin condition.

FDA cleared for the non-invasive treatment of wrinkles and rhytids, including the upper and lower eyelids. FDA cleared for the temporary improvement in the appearance of cellulite when combined with vibration.

Side effects are generally mild and may include redness, swelling, blisters, bumps and rare instances of surface irregularities.

© 2009 Solta Medical, Inc. All rights reserved. Thermage, Thermage CPT and Comfort Pulse Technology are trademarks or registered trademarks of Solta Medical, Inc.

MK2086